



MEDICAL INDUSTRIES

Hoffmann Family of Companies



Elite II™

NOA's Elite II low bed offers your facility a variety of standard features, all at an affordable price.

Standard Features

- Travels from a low of 7" to a convenient high height of 30"
- Center step safety lock
- 600 lb. safe working load
- Self-leveling motors
- 4 Year warranty on electronics and 15 year warranty on bed deck and frame
- Swivel locking casters and flip down caster guides
- Baseboard bumper
- 10 Button handset
- Meets 60601-2-52

Optional Features

- Head and footboards- light oak, wild cherry, walnut, mahogany
- Soft touch side rails
- Standard or folding assist bar
- Pressure redistribution foam mattress
- Footboard mounted staff control
- Trapeze



MEDICAL INDUSTRIES
Hoffmann Family of Companies



Technical Specifications

Elite II™

Bed length	80"
Bed width	35"
Bed width, with side rails	40"
Overall bed length with head & foot boards	84.5"
Deck height range	7" - 30"
Safe working load	600 lbs.
Maximum back angle	65°
Maximum hip-knee angle	30°



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com