



MEDICAL INDUSTRIES

Hoffmann Family of Companies



## NOA Light™

The NOA Light offers a quality, cost effective high low bed that meets the needs of both the resident and clinical staff.

### Standard Features

- Rises from a low of 8.5" to a high of 25"
- 4 DC Motors that provide elevation, head and foot adjustment
- Sturdy slat deck providing a solid sleep surface and mattress ventilation
- Standard sleep deck is 35" wide by 80" long
- Four locking casters with two directional casters
- Move in any position
- Three year warranty on electronics and 10 year warranty on bed, deck, and frame

### Optional Features

- Safety lock for immobilization
- Head and footboards- light oak, wild cherry, walnut, mahogany
- Soft touch side rails
- Standard or folding assist bar
- Pressure redistribution foam mattress
- Footboard mounted staff control
- Trapeze
- Deck expansion kit



## Technical Specifications

### NOA Light

Bed length .....	80"
Bed width .....	35"
Bed width, with side rails .....	40"
Overall bed length with head & foot boards .....	84.5"
Deck height range .....	8.5" - 25"
Safe working load .....	450 lbs.
Maximum back angle .....	60°
Maximum hip-knee angle .....	40°

