



Platinum ML Hospital Bed

Economically Priced with Top of the Line Features

Features

- Manual CPR quick release head deck
- Under bed lighting
- Battery backup
- High impact quick release head and footboards
- High impact side rails
- Central and directional locking system with 5" casters
- Transport handles for ease of transporting bed and patient
- 15 year warranty on deck & frame welds
- 4 year warranty on electronics
- Drainage hooks located on both sides of bed frame
- Caster base covers
- Self leveling motors
- 600 lb. safe working load
- Foot board control



MEDICAL INDUSTRIES

Hoffmann Family of Companies



Technical Specifications

Platinum ML Hospital Bed

| | | | |
|---------------------------------------|------------------------|--------------------------------------|-----------|
| Frame width with side rails | 42" | Bed clearance | 5.5" |
| Frame length with head and footboard | | Safe working load | 600 lbs. |
| Standard length | 84.5" | Deck height range | 11" - 31" |
| Sleep deck | | Max & reverse trendelenburg. | 14° |
| Width | 35" | Max back angle | 68° |
| Length | Available in 80" - 84" | Max hip-knee angle | 45° |



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com