



**MEDICAL INDUSTRIES**

Hoffmann Family of Companies



## EliteCare

NOA Medical's EliteCare is our most cost-effective bed to date. Fully functional with a three-motor system. Designed for your convenience.

### Standard Features

- Three motor system
- Travels from a low of 7.75" To a high of 29"
- Self-leveling
- 8 Button handset
- 450 lbs. safe working load
- Swivel locking casters
- Removable baseboard bumper
- Grid deck design
- 3 Year electronic warranty & 10 year bed, deck, frame warranty
- Meets 60601-2-52

### Optional Features

- Head & footboards - light oak, wild cherry, walnut, mahogany
- Soft touch assist rail
- Standard or folding assist bar
- Pressure redistribution foam mattress
- Trapeze



MEDICAL INDUSTRIES

Hoffmann Family of Companies



## Technical Specifications

### Elite Care

Bed length .....	80"
Bed width .....	35.5"
Overall bed length with head & foot boards .....	85"
Deck height range .....	7.75" - 29"
Max back angle .....	72°
Max hip-knee angle .....	43°
Safe working load .....	450 lbs.



801 Terry Lane, Washington, MO 63090  
800.633.6068 • noamedical.com