



## NOA Platinum Patient Transport Stretcher

Our hydraulic patient transport stretcher offering you multiple standard features at an affordable price.

### Features

- Powder coated steel frame and ABS base
- Back board is X-Ray translucent and is fitted with X-Ray cassette holder
- Foldable 5-staff aluminum alloy guardrail
- Equipped with two USA hydraulic pumps
- Five functions: adjustable back rest, leg rest, trendelenburg/reverse trendelenburg and height adjustment
- Telescoping IV pole
- 5th wheel assist for directional control during transport
- 8 inch casters with central locking
- Two retractable push handles for controlled movement
- Mattress with moisture resistant cover
- Oxygen cylinder holder/tray for placement in ABS base
- Rubber corner bumpers
- 1 year warranty



## Technical Specifications

### **NOA Platinum Patient Transport Stretcher**

Overall stretcher width . . . . .	37.5"
Overall stretcher length . . . . .	84"
Safe working load . . . . .	600 lbs.
Deck height range . . . . .	24.5" - 36"
Max & reverse trendelenburg . . . . .	-10° / +10°
Caster size . . . . .	8"
Max back angle . . . . .	70°
Max hip-knee angle . . . . .	28°

