


Protekt® 500

Gel Infused Foam Pressure Redistribution Mattress



FEATURES AND BENEFITS:

The Protekt® 500 with (VGPT)™ Visco Gel Polymer Technology is a revolutionary support surface that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam.

- Reduces occurrence of pressure sores by channeling heat away from the body significantly faster than traditional foam, virtually eliminating perspiration.
- Superior pressure redistribution is delivered by providing more support and less cradling at higher levels of compression, redistributing pressure over a large surface area eliminating static fatigue.
- Protects vulnerable heel area with a 30 degree slope redistributing pressure to the entire leg and away from the heel.
- PTSS Rail System provides additional support for easy patient transfer and edge of bed sitting.
- Stretch nylon top cover reduces perspiration and

provides a low shear friction surface that is anti-microbial and fluid resistant.

- Bottom Cover: Triple laminate vinyl; fluid resistant; anti-microbial/anti-bacterial; fire retardant.
- Available in 42", 48", 54" & 60" widths.
- Fire sleeve meets Federal Fire Code 16 CFR part 1633.
- Weight Capacity: 450 lbs.
- Warranty: Core; 15 Year Non-Prorated. Cover; 2 Year.
- Optional raised rails available on all sizes.

81051	Protekt® 500	36"x76"x6"
81052	Protekt® 500	36"x80"x6"
81053	Protekt® 500	36"x84"x6"
81054	Protekt® 500 Plus	36"x76"x9"
81055	Protekt® 500 Plus	36"x80"x9"
81056	Protekt® 500 Plus	36"x84"x9"



801 Terry Lane, Washington, MO 63090

800.633.6068 www.noamedical.com