



## Elite Riser™

### Four Year Warranty on Electronics

Lifetime warranty on bed deck and frame.

### 600 lb. safe working load

### Faster & Quieter

The Elite Riser swiftly moves from 7.5" to 28".

### Lightweight

The Elite Riser is lighter than most high/low beds.

### Self Leveling

A built-in sensor in the control box communicates with the actuators to keep the bed level at all times.

### Very Mobile

The Elite Riser can be easily moved at any height and in any position.

### Safety Lock Option

Makes it one of the safest beds on the market.

# NOA Elite Riser™

## Features Include:

- Bed goes straight up and straight down
- Swivel braking casters at all four corners
- Bed rolls and locks at any height
- Entrapment compliant with properly sized mattress
- Easy to maintain
- Faster
- Quieter

## Optional Features:

- Safety lock for immobilization
- Head and Foot Boards - Light Oak, Dark Oak, Wild Cherry, Walnut, Mahogany
- Soft-Touch Side Rails
- Sub-acute handset
- Assist Bar, Left, Right and Universal
- Floor Pad, 2" thick, x 2' x 6'
- Therapeutic Foam Mattresses
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- High Impact Head and Foot Boards

ITEM	76" length	80" length	84" length
bed width	35"	35"	35"
bed width, with side rails	40.0"	40.0"	40.0"
bed width, with assist bar	36.5"	36.5"	36.5"
overall bed length with head and foot boards	80.5"	84.5"	88.5"
deck height range	7.5" to 28"	7.5" to 28"	7.5" to 28"
safe working load	600 lbs	600 lbs	600 lbs



**MEDICAL INDUSTRIES INC.**

801 Terry Lane Washington, MO 63090  
 Toll Free: 1-800-633-6068 | [www.noamedical.com](http://www.noamedical.com)