



NOA R600

Self-Leveling

A built-in sensor in the control box communicates with the actuators to keep the bed level at all times.

600 lb. Safe Working Load

4 Year Warranty

on electronics and lifetime warranty on bed and deck frame.

Hospital Grade Casters

allow for easy transport within the facility

Fully Adjustable

The R600 provides trendelenberg positioning as well as cardiac chair positioning.

The R600 Bed Comes Complete

with soft touch side rails and impact resistant polyethylene head and foot boards making it a perfect addition to your LTC facility or rental fleet of bed systems.

NOA R600

Features Include:

- Height Range from 9.5" low to 29.5" high
- Hospital Grade 4" Locking Casters
- High density, white Polyethylene head and foot boards
- Rib Deck design for added durability and increased mattress ventilation
- Can be moved in any position and at any height
- Goes straight up and straight down

Optional Features:

- Soft Touch Side Rails
- Footboard Mounted Staff Control
- In-Rail Pendant
- Assist Bar, Left, Right and Universal
- Central Locking System (locks all four casters with one motion)
- Lock Out Hand Set
- High Impact Plastic head and foot boards- Standard Laminate Colors Available
- Full Overbed Trapeze
- Floor Pad, 2" thick, x 2'x6'

ITEM	76" length	80" length	84" length
bed width	35"	35"	35"
bed width, with side rails	40.0"	40.0"	40.0"
bed width, with assist bars	36.5"	36.5"	36.5"
overall bed length			
with head and foot boards	80.5"	84.5"	88.5"
deck height range	9.5" to 29.5"	9.5" to 29.5"	9.5" to 29.5"
safe working load	600 lbs	600 lbs	600 lbs



NOA

MEDICAL INDUSTRIES INC.

801 Terry Lane Washington, MO 63090
 Toll Free: 1-800-633-6068 | www.noamedical.com