



NOA Twin Riser

The NOA Twin Riser Bed incorporates an extra wide sleep surface for a new level of resident comfort and relaxation that makes them feel like they're back in their own bed. No more tossing and turning in standard size long-term care beds.

Wide Sleep Deck

All NOA Twin Riser beds come with an extra-wide 39" sleep surface

Touch a Button

And watch the wide NOA Twin Riser gently rise from 7.5" to 28" for resident safety and caretaker convenience.

Touch Another Button

And watch as the Twin Riser adjusts the head and leg sections of the bed to fit your residents needs, and/or easily adjusts to a comfort chair position.

Guarantee a Good Nights Sleep

With 39" of sleep space your residents will have more freedom to move and turn in bed.

NOA Twin Riser

Features Include:

- A 39" wide sleep surface for superior comfort
- High-Low function from 7.5" to 28" for resident safety and caregiver convenience
- 600 pound safe working load
- Four quiet DC motors for smooth operation
- Effortless mobility in high or low positions
- Long linen life, excellent mattress aeration, and easy cleaning
- Four Year Warranty on Electronics
Lifetime warranty on bed deck and frame

Optional Features:

- Safety lock for immobilization
- Head and Foot Boards- Light Oak, Dark Oak, Wild Cherry, Walnut, Mahogany
- Soft-Touch Side Rails
- Sub-acute handset
- Assist Bar, Left, Right and Universal
- Floor Pad, 2"thick, x 2'x6'
- Therapeutic Foam Mattress
- Footboard Mounted Staff Control
- Full Overbed Trapeze

| ITEM | 80" length | 84" length |
|-----------------------------|-------------|-------------|
| bed width | 39" | 39" |
| bed width, with side rails | 44.0" | 44.0" |
| bed width, with assist bars | 40.5" | 40.5" |
| overall bed length | | |
| with head and foot boards | 84.5" | 88.5" |
| deck height range | 7.5" to 28" | 7.5" to 28" |
| safe working load | 600 lbs | 600 lbs |



NOA

MEDICAL INDUSTRIES INC.

801 Terry Lane Washington, MO 63090

Toll Free: 1-800-633-6068 | www.noamedical.com