



NOA Light™

NOA Medicals newest high/low bed system quickly rises from a low of 8.5" to a high of 25".

Fully adjustable

The four motors provides your residents with comfort positioning at the head and foot area, as well as height adjustment.

Head & Footboards

Available with standard wood grain head & footboards or High-Impact head and footboards with quick release

Effectively meets your needs

The NOA Light meets the needs of your clinical staff as well. The NOA Light provides a safe working load of 450 lbs. and offers a number of quality standard accessories as well.

Slat Deck Design

The NOA Light offers a new slat deck design to provide a high quality, solid sleep surface for your resident.

NOA Light™

Features Include:

- Rises from a low of 8.5" to high of 25"
- Has 4 DC Motors that provide elevation, head and foot adjustment.
- Has a sturdy slat deck providing a solid sleep surface and mattress ventilation
- 35" wide by 80" long
- Four Locking Casters with two directional casters
- Can be moved in any position
- Easy to clean and disinfect
- Three year warranty on electrical and mechanical parts. Fifteen year warranty on bed deck and frame.

Optional Features:

- Soft Touch Side Rails
- Assist Bar, Left, Right, Universal
- Footboard Mounted Staff Control
- Standard wood grain head and foot boards
- Floor Pad, 2" thick, x 2'x6'
- Safety Locks
- Expansion Kit (expands 35" deck to 39" & 42")

ITEM

80" length

bed width	35"
bed width, with side rails	40.0"
bed width, with assist bars	36.5"
overall bed length	
with head and foot boards	84.5"
deck height range	8.5" to 25"
safe working load	450 lbs



NOA

MEDICAL INDUSTRIES INC.

801 Terry Lane | Washington, MO 63090
Toll Free: 1-800-633-6068 | www.noamedical.com