



## NOA Bariatric Bed

### Wide Sleep Deck

All NOA Bariatric Beds come with a comfortable 42" wide sleep deck that easily adjusts from 80" to 84" in length.

### 800 lb Safe Working Load

Addresses a majority of your bariatric resident population weight requirements. 800 lb. safe working load allows even your heaviest residents to benefit from the variable heights and positions standard with the NOA Bariatric Bed.

### Lock it In Place

The Bariatric Bed can be moved at any height, at any position or locked in place.

### Recliner Position

For added comfort, the bed can be put into a Recliner Position with the bed in a Reverse Trendelenberg position and the knee gatch and head deck raised for greater resident comfort.

### Trolley System

NOA Bariatric Beds come with a trolley system which allows movement of the bed from one room to the other or upright storage to save space when not in use.

### New Options

- In Bed Scale with digital readouts
- Nurse Control Panel
- CPR Quick Release

# NOA Bariatric Bed

## Features Include:


- Soft-Touch Side Rails
- Polyethylene Head and Foot Boards
- Transport Trolley System
- T/RT Capable
- Baseboard Bumper
- Entrapment Compliant
- Four year warranty on electronics.  
Lifetime warranty on bed deck and frame.

## Optional Features:

- Lockout Handset
- Universal Assist Bar
- Footboard Mounted Staff Control
- Therapeutic Foam Mattress
- Full Overbed Trapeze
- Floor Pad, 2" thick, x 2'x6'
- In Bed Scale
- CPR Quick Release
- Out of Bed Alarm & Nurse Call
- Under Bed Light
- Battery Back Up

ITEM	80" length	84" length
bed width	42.0"	42.0"
bed width, with side rails	47.0"	47.0"
overall bed length		
with head and foot boards	84.0"	88.5"
deck height range	9.5" to 30"	9.5" to 30"
safe working load	800 lbs	800 lbs

Trolley System



# NOA

MEDICAL INDUSTRIES INC.

801 Terry Lane | Washington, MO 63090  
Toll Free: 1-800-633-6068 | [www.noamedical.com](http://www.noamedical.com)