



Elite EX™

Part of the Elite Series of Beds

Integrated Deck Expansion System

Easily expand your sleep surface from 36" to 39" to 42" width through the simple pull of a button, and 80" to 84" length.

Four Year Warranty on Electronics

15 year warranty on bed deck and frame welds.

Low Height

The Elite EX travels from a low of 7" to a convenient high height of 30".

Center Step Safety Lock Comes Standard

Lock your bed in place with one quick step.

Fully Adjustable

The four motors provide your residents with comfort positioning at the head and foot area, as well as height adjustment.

600 lbs. Safe Working Load

Elite EX™

Part of the Elite Series of Beds

Features Include:

- Integrated deck expansion system
- 12 button handset
- Trapeze holder located at the head of bed
- Center step safety lock
- Wall bumpers
- Swivel locking casters
- Flip down caster guides for ease in steering

Optional Features:

- Standard wood grain head and foot boards
- Soft touch side rails
- Universal or folding assist bar
- Foot board mounted staff control
- Trapeze
- Therapeutic foam mattress
- Floor pad, 2" thick, x 2' x 6'

ITEM

Bed Width

Bed width, with side rails:

Overall bed length with
head and foot boards:

Deck height range:

80" Length (expands easily to 84")

35" (can be expanded to 39" or 42")

40"

84.5"

7"-30"



NOA

MEDICAL INDUSTRIES INC.

801 Terry Lane Washington, MO 63090
Toll Free: 1-800-633-6068 | www.noamedical.com