



Elite II

Part of the Elite Series of Beds

Center Step Safety Lock Comes Standard
Lock your bed in place with one quick step.

Four Year Warranty On Electronics
15 year warranty on bed deck and frame welds.

Low Height

The **Elite II** travels from a low of 7" to a convenient high height of 30".

Fully Adjustable

The four motors provide your residents with comfort positioning at the head and foot area, as well as height adjustment.

600 lbs. Safe Working Load

Convenient 12 Button Handset
Auto Contour & Comfort Chair feature

Elite II

Part of the Elite Series of Beds

Features Include:

- 12 button handset
- Center step safety lock
- Swivel locking casters
- Flip down caster guides for ease in steering


Optional Features:

- Standard wood grain head and foot boards
- Soft touch side rails
- Universal or folding assist bar
- Foot board mounted staff control
- Therapeutic foam mattress
- Beveled floor pad, 2" thick x 2' x 6'

Item

80" Length

Bed width	35"
Bed width, with side rails	40"
Overall bed width with head and foot boards	84.5"
Deck height range	7"-30"



801 Terry Lane Washington, MO 63090
Toll Free: 1-800-633-6068 www.noamedical.com