



Elite II™

NOA's Elite II low bed offers your facility a variety of standard features, all at an affordable price.

Standard Features

- Travels from a low of 7" to a convenient high height of 30".
- Unique anti-shear head deck, easily keeps resident from sliding down in bed.
- Center Step Safety Lock
- 600 lb. safe working load
- Self- Leveling motors
- Four year warranty on electronics and 15 year warranty on bed deck and frame
- Swivel locking casters and flip down caster guides
- Baseboard bumper
- 12 button handset

Optional Features

- Head and Footboards- Light Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Standard or Folding Assist Bar
- Pressure Redistribution Foam Mattress
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- Fall Mat



Technical Specifications

Bed length	80"
Bed width	35"
Bed width, with side rails	40"
Overall bed length with head & foot boards	84.5"
Deck height range	7" - 30"