



NOA Light™

The NOA Light offers a quality, cost effective high low bed that meets the needs of both the resident and clinical staff.

Standard Features

- Rises from a low of 8.5" to a high of 25"
- 4 DC Motors that provide elevation, head and foot adjustment
- Sturdy slat deck providing a solid sleep surface and mattress ventilation
- Standard sleep deck is 35" wide by 80" long
- Four locking casters with two directional casters
- Move in any position
- Three year warranty on electrical and mechanical parts.
- Fifteen year warranty on bed deck and frame.

Optional Features

- Safety lock for immobilization
- Head and Footboards- Light Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Standard or Folding Assist Bar
- Pressure Redistribution Foam Mattress
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- Deck Expansion Kit
- Fall Mat



Technical Specifications

Bed length	80"
Bed width	35"
Bed width, with side rails	40"
Bed width, with assist bar	36.5"
Overall bed length with head & foot boards.....	84.5"
Deck height range	8.5" - 25"
Safe working load	450 lbs.