



Twin Riser™

Incorporates an extra wide sleep surface for a new level of resident comfort and relaxation that makes them feel like they're back in their own bed. No more tossing and turning in standard size long-term care beds.

Standard Features

- A 39" wide sleep surface for superior comfort
- High-Low function from 7.5" to 28"
- Effortless mobility and moves straight up and down at any height and in any position
- 600 lb. safe working load
- Self-Leveling motors
- Four year warranty on electronics and lifetime warranty on bed deck and frame
- Swivel braking casters at all four corners
- Entrapment compliant
- Slat deck design provides long linen life, excellent mattress aeration and easy cleaning

Optional Features

- Safety lock for immobilization
- Head and Footboards- Light Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Standard or Folding Assist Bar
- Sub-Acute handset
- Pressure Redistribution Foam Mattress
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- Deck Expansion Kit
- Fall Mat



Technical Specifications

Bed length	80"	84"
Bed width	39"	39"
Bed width, with side rails	44"	44"
Bed width, with assist bar	40.5"	40.5"
Overall bed length with head & foot boards.....	84.5"	88.5"
Deck height range	7.5" - 28"	7.5" - 28"
Safe working load	600 lbs.	600 lbs.