

## What do residents everywhere say about this **extra-wide** bed?



**The NOA Twin Riser™ Bed incorporates an extra-wide sleep surface for a new level of resident comfort and relaxation that makes them feel like they're back in their own twin bed. No more tossing and turning in standard-size resident beds.**

#### **The NOA Twin Riser delivers:**

- A 39" wide sleep surface for superior comfort
- High-Low function from 7-26" for resident safety and caregiver convenience
- 500 pound safe working load for larger residents
- Four quiet DC motors for smooth operation
- Effortless mobility in up or down positions
- Long linen life, excellent mattress aeration and easy cleaning

Call **800-633-6068** today.

[www.noamedical.com](http://www.noamedical.com)



# TWIN RISER