


Geo-Mattress® 350

6" Therapeutic Foam Mattress



An effective solution
for pressure redistribution
in all care settings.

Item*	Description
6000008	Geo-Mattress 350, 75"L x 36"W x 6"H
6000009	Geo-Mattress 350, 80"L x 36"W x 6"H
6000010	Geo-Mattress 350, 84"L x 36"W x 6"H
6000012	Geo-Mattress 350, 80"L x 39"W x 6"H
6000013	Geo-Mattress 350, 84"L x 39"W x 6"H

Weight Limit: 350 lbs.

Warranty: 5 years, (2 years pro-rated).

Flammability Protection:
includes inner fire barrier sleeve
for certified compliance to NFPA 101
(Life Safety Code™); Cal Tech Bulletin #117,
16 CFR 1632 and 1633; and BFD IX-11.

Cross-cut top surface
creates 800 individual
cells in 3 zones for en-
hanced immersion and
pressure redistribution.



Fluid-proof, wipe-clean
coated nylon cover, with
hidden bottom zipper for
easier, more thorough
cleaning.



Manufactured for:



WWW.NOAMEDICAL.COM 800.633.6068

SPANAmerica
Innovative Solutions