



## Elite EX™

NOA Medical's Elite EX leads the long term care industry in providing a wider, more comfortable sleep deck for residents.

### Standard Features

- Integrated deck expansion system, toolless expansion of bed from 35" to 39" to 42" width through a simple pull of a knob, and 80" to 84" in length
- Travels from a low of 7" to a convenient high height of 30"
- Self-leveling motors
- 10 button handset
- 600 lb. safe working load
- Center step safety lock
- Swivel locking casters and flip down caster guides
- Baseboard bumper
- 4 Year warranty on electronics and 15 year warranty on bed deck and frame
- Meets 60601-2-52

### Optional Features

- Head and footboards - light oak, wild cherry, walnut, mahogany
- Soft touch assist rail
- Standard or folding assist bar
- Pressure redistribution foam mattress
- Footboard mounted staff control
- Trapeze



MEDICAL INDUSTRIES

Hoffmann Family of Companies



## Technical Specifications

### Elite EX™

Bed length (Extends w/o tools) . . . . .	80" - 84"
Bed width (Extends w/o tools) . . . . .	35" - 39" - 42"
Overall bed length with head and foot boards . . . . .	84.5" - 88.5"
Deck height range . . . . .	7" - 30"
Maximum back angle . . . . .	65°
Maximum hip-knee angle . . . . .	30°
Safe working load . . . . .	600 lbs.



801 Terry Lane, Washington, MO 63090  
800.633.6068 • noamedical.com