

Elite Riser Bed swiftly moves straight up and straight down from 7.5" to 28", powered by whisper quiet motors.

## **Standard Features**

- Moves straight up and down, and can be moved at any height and in any position
- · 600 lb. safe working load
- · Self-Leveling motors
- Four year warranty on electronics and lifetime warranty on bed deck and frame
- · Swivel braking casters at all four corners
- · Entrapment compliant
- · Fast and Quiet
- · Easy to maintain

## **Optional Features**

- · Safety lock for immobilization
- Head and Footboards- Light Oak, Wild Cherry, Walnut, Mahogany
- · Soft Touch Side Rails
- · Standard or Folding Assist Bar
- · Sub-Acute handset
- · Pressure Redistribution Foam Mattress
- · Footboard Mounted Staff Control
- Full Overbed Trapeze
- · Deck Expansion Kit
- · Fall Mat



## **Technical Specifications**

Bed length	76"	80"	84"
Bed width	35"	35"	35"
Bed width, with side rails	40"	40"	40"
Bed width, with assist bar	36.5"	36.5"	36.5"
Overall bed length with			
head & foot boards	80.5"	84.5"	88.5"
Deck height range	7.5" - 28"	7.5" - 28"	7.5" - 28"
Safe working load	600 lbs.	600 lbs.	600 lbs.

