

Overbed Tables

NOA Medical, a leader in the development and manufacture of bed systems and quality accessories for the patient/resident room, introduces our newly designed Over Bed Table. Our new OBT provides increased stability and durability to your residents through a unique U-Shaped base while providing unobstructed movement around day chairs and resident beds.

New & Improved:

- High pressure laminate top measures 15" X 31.5"
- Table tops come in five standard wood grains, i.e. Light Oak, Dark Oak, Walnut, Mahogany, and Wild Cherry. Custom colors can be quoted upon request.
- Adjustable height ranging from a low of 29" to a high of 43" for the added convenience of your resident.
- Base is powdercoated to aid in cleaning and infection control.


New Base Design:

The new base uses rectangular tubing for added strength, adding to the durability and stability of the table. Rated for 50 lb weight capacity.



Other Quality Accessories from NOA Medical

Assist Bar



Mattresses



Deck Expansion System

