



Platinum NS Hospital Bed

Meeting the needs of your patients and
Clinical Staff at an Affordable Price

Features

- One Button CPR quick release head deck
- Under bed lighting
- Battery backup
- High impact quick release head and footboards
- High impact side rails
- Central and directional locking system with 5" Casters
- Transport handles for ease of transporting bed and patient
- 15 Year warranty on deck & frame welds
- 4 Year warranty on electronics
- Drainage hooks located on both sides of bed frame
- Self leveling motors
- 600 Lb. Safe working load
- Digital nurse panel with real time readouts
- Meets 60601-2-52



Technical Specifications

Platinum NS Hospital Bed

Frame width with side rails	42"	Bed clearance	5.5"
Frame length with head and footboard		Safe working load	600 lbs.
Standard length	84.5"	Deck height range	11" - 31"
Sleep deck		Max & reverse trendelenburg . . .	14°
Width	35"	Max back angle	70°
Length	Available in 80" - 84"	Max hip-knee angle	35°

