



## Overbed Tables



### Standard Overbed Table Features

- New thermofoil non drip top measures 15" x 31.75"
- Tabletops comes in four standard wood grains; light oak, walnut, mahogany and wild cherry. Custom colors available
- Adjustable height ranging from a low of 29" to a high of 43" for the added convenience of your resident
- The U-base design uses rectangular tubing for added strength, adding to the durability stability of the table
- Base is available in black or beige and is rated for a 50 lb. weight capacity
- Base is powder coated to aid in cleaning and infection control
- Thermofoil top has a five year warranty
- Base and mechanical parts have a one year warranty

### Gas Cylinder Overbed Table Features

- High pressure laminate top measures 15" x 31.75"
- Table tops comes in two standard wood grains; light oak or mahogany
- Adjustable height ranging from a low of 29" to a high of 43" for the added convenience of your resident
- The U-Base design uses rectangular tubing for added strength, adding to the durability stability of the table
- Base is available in beige and is rated for a 50 lb. weight capacity
- Base is powder coated to aid in cleaning and infection control
- Top, base and mechanical parts have a one year warranty