



MEDICAL INDUSTRIES

Hoffmann Family of Companies



NOA Bariatric Elite Bed

The NOA Bariatric Elite Bed will accommodate virtually all of the bariatric bed requirements in your facility.

Standard Features

- 800 lb. Safe Working Load
- High-low Function From 9.5"-30"
- 10 Button Handset Includes T/RT
- Swivel Locking Casters
- Bed Rolls and Locks at any Height
- Mattress Retainers
- 4 Year Warranty on Electronics
- 15 Year Warranty on Bed, Deck, and Frame
- Meets 60601-2-52

Optional Features

- HDPE Head & Foot Boards
- Soft Touch Side Rails
- Standard & Folding Assist Bars
- Underbed Lighting
- Drainage Hooks
- CPR Quick Release
- USB Ports
- Footboard Mounted Staff Control
- Therapeutic Foam Mattress
- Fall Mat
- Freestanding Trapeze
- Digital Nurse Panel
- Transport Trolley System



MEDICAL INDUSTRIES
Hoffmann Family of Companies



Technical Specifications

NOA Elite Bariatric Bed

	80" Length	84" Length
Bed width, with side rails	46"	46"
Sleep deck length	80"	84"
Overall bed length	84.5"	88.5"
Deck height range	9.5" to 30"	9.5" to 30"
Max back angle	65°	65°
Max hip-knee angle	45°	45°
Safe working load	800 lbs.	800 lbs.
Length expansion kit included		



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com