



NOA Expandable Bariatric Hospital Bed with In-Bed Scale

The NOA Expandable Bariatric Hospital Bed addresses all your bariatric needs allowing even your heaviest residents variable heights and positions as standard.

Standard Features

- Toolless Expandable Bed Width - 42" - 45" - 48"
- 1000 lb. Safe Working Load
- In Bed Scale System
- Out of Bed Alarm that Ties into the Nurse Call System
- T/RT Capability
- Entrapment Compliant
- Central & Directional Locking System
- 4 Year Warranty on Electronics
- 15 Year Warranty on Bed Deck & Frame
- Meets 60601-2-52

Optional Features

- HDPE Head & Foot Boards
- Soft Touch Side Rails
- Standard & Folding Assist Bars
- Therapeutic Foam Mattress
- Fall Mat
- Battery Backup
- USB Ports
- Underbed Lighting
- CPR Quick Release
- Digital Nurse Panel
- Freestanding Trapeze
- Transport Trolley System



MEDICAL INDUSTRIES

Hoffmann Family of Companies



Technical Specifications

NOA Expandable Deck Hospital Bariatric Bed

Bed length	80"
Bed width	42" - 45" - 48"
Overall bed length with head & foot boards	84"
Deck height range	11.5" - 31.75"
Max back angle	65°
Max hip-knee angle	45°
Safe working load	1000 lbs.
Length expansion kit included	



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com