



## NOAH Hospital Bariatric Bed with In Bed Scale

The NEW NOAH Bariatric Hospital Bed addresses all your bariatric needs allowing even your heaviest residents variable heights and positions as standard.

- High Impact Side Rails
  - HDPE Head and Foot boards
  - Transport trolley system
  - 1000 lb. safe working load
  - T/RT Capability
  - Entrapment compliant
  - Four year warranty on electronics.
  - Lifetime warranty on bed deck and frame.
- Central and Directional Locking Casters
  - Digital Nurse Panel
  - Underbed Light
  - CPR Quick Release
  - In Bed Scale System
  - Out of Bed Alarm that ties into the nurse call system
- Optional Features:**

  - Universal or Folding Assist Bar
  - Therapeutic Foam Mattress
  - Full Overbed Trapeze
  - Floor Pad, 2" thick, x 2'x6'
  - Battery Back Up



## Technical Specifications

Bed length	80"	84'''
Bed width . . . . .	42"	42"
Bed width, with side rails . . . . .	47"	47"
Overall bed length with head and foot boards. . . . .	84"	88.5"
Safe working load . . . . .	1000 lbs.	1000 lbs.
Deck height range: . . . . .	11.50" to 31.75"	11.50" to 31.75"
Optional Deck Expansion System Bed width. . . . .	48"	48"