



MEDICAL INDUSTRIES

Hoffmann Family of Companies



NOA Expandable Deck Hospital Bariatric Bed with In Bed Scale

The NOA Bariatric Hospital Bed addresses all your bariatric needs allowing even your heaviest residents variable heights and positions as standard.

Standard Features

- Toolless expandable bed width - 42" - 45" - 48"
- High impact side rails
- HDPE head & foot boards
- In bed scale system
- Out of bed alarm that ties into the nurse call system
- 1000 lb. safe working load
- T/RT capability
- Entrapment compliant
- Four year warranty on electronics
- Lifetime warranty on bed deck & frame
- Central & directional locking system
- Digital nurse panel
- Underbed light
- CPR quick release

Optional Features

- Universal or folding assist bar
- Therapeutic foam mattress
- Full overbed trapeze
- Floor pad, 2" thick, x 2' x 6'
- Battery back up



MEDICAL INDUSTRIES

Hoffmann Family of Companies



Technical Specifications

NOA Expandable Deck Hospital Bariatric Bed

- Bed length 80"
- Bed width 42" - 45" - 48"
- Overall bed length with head & foot boards 84"
- Deck height range 11.5" - 31.75"
- Max back angle 65°
- Max hip-knee angle 45°
- Safe working load 1000 lbs.
- Length expansion kit included



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com