



## NOA Bariatric Hospital Bed

The NOA Bariatric Hospital Bed addresses all your bariatric needs allowing even your heaviest residents variable heights and positions as standard.

### Standard Features

- Toolless expandable deck bed width - 42" - 45" - 48"
- High impact side rails
- HDPE head & foot boards
- 1000 lb. safe working load
- T/RT capability
- Entrapment compliant
- Four year warranty on electronics
- Lifetime warranty on bed deck & frame
- Central & directional locking system
- Digital nurse panel
- Underbed light
- CPR quick release

### Optional Features

- Universal or folding assist bar
- Therapeutic foam mattress
- Full overbed trapeze
- Floor pad, 2" thick, x 2' x 6'
- Battery back up



Trolley System



## Technical Specifications

### NOA Hospital Bariatric Bed

Bed length . . . . .	80"
Bed width . . . . .	42"
Overall bed length with head & foot boards . . . . .	84"
Deck height range . . . . .	11.5" - 31.75"
Expandable deck . . . . .	42" - 45" - 48"
Max back angle . . . . .	65°
Max hip-knee angle . . . . .	45°
Safe working load . . . . .	1000 lbs.
4" length extension kit included	

