



EliteCare

NOA Medical's EliteCare is our most cost-effective bed to date. Fully functional with a three-motor system. Designed for your convenience.

Standard Features

- Three motor system
- Travels from a low of 7.75" To a high of 29"
- Self-leveling
- 8 Button handset
- 450 lbs. safe working load
- Swivel locking casters
- Removable baseboard bumper
- Grid deck design
- 3 Year electronic warranty & 10 year bed, deck, frame warranty
- Meets 60601-2-52
- Trapeze

Optional Features

- Head & footboards - light oak, wild cherry, walnut, mahogany
- Soft touch assist rail
- Standard or folding assist bar
- Pressure redistribution foam mattress
- Fall mat



Technical Specifications

Elite Care

| | |
|--|-------------|
| Bed length | 80" |
| Bed width | 35.5" |
| Overall bed length with head & foot boards | 85" |
| Deck height range | 7.75" - 29" |
| Max back angle | 72° |
| Max hip-knee angle | 43° |
| Safe working load | 450 lbs. |

