

NDA

MEDICAL INDUSTRIES INC.

ELITE FAMILY OF PRODUCTS





About Us

For over 30 years NOA Medical Industries, Inc. has been a recognized leader in the development and manufacture of high quality, fully adjustable, electric and manual high / low bed systems and sleep safety accessories for the long-term care market. We pride ourselves on our ability to continually develop new bed systems and sleep safety products that effectively meet the needs of our customers and their residents. NOA Medical has quickly become the standard from which to measure high / low bed systems, their quality, function and design. And NOA Medical is not stopping there. From our combined 90,000 square foot facilities in Washington, MO, NOA Medical continues to design and innovate, and has plans to introduce a number of new, market driven bed systems and accessories over the coming year.

NOA Medical. Innovators in Quality, Safety and Design.

The
NOA *Elite* Family
Only the best can be called Elite



Table of Contents

About Us	3
Beds	
Elite Riser	6
Twin Elite Riser	8
Elite Bariatric	10
Elite R600	12
LTC Light	14
Accessories	
Standard Head and Foot Boards	16
Impact Resistant Head and Foot Boards	16
Side Rails	16
Sub Acute Handset	16
In Rail Pendant	17
Safety Lock	17
Assist Bar	17
Floor Mats	17
Over Bed Table	18
Deck Expansion System	18
Trapeze	19

NOA Elite Riser™

Four Year Warranty on Electronics
Lifetime warranty on bed deck and frame.

600 lb. safe working load

Faster & Quieter

The new Elite Riser swiftly moves from 7.5" to 28".

Lightweight

The Elite Riser is lighter than most high/low beds.

Self Leveling

A built-in sensor in the control box communicates with the actuators to keep the bed level at all times.

Very Mobile

The Elite Riser can be easily moved at any height and in any position.

New Safety Lock Option

Makes it one of the safest beds on the market.



Features Include:

- Bed goes straight up and straight down
- Swivel braking casters at all four corners
- Bed rolls and locks at any height
- Entrapment compliant with properly sized mattress
- Easy to maintain
- Faster
- Quieter

Optional Features:

- Safety lock for immobilization
- Head and Foot Boards - Light Oak, Dark Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Sub-acute handset
- Assist Bar, Left, Right and Universal
- Floor Pad, 2" thick, x 2' x 6'
- Therapeutic Foam Mattresses
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- High Impact Plastic head and foot boards
Standard Laminate Colors available

ITEM	76" length	80" length	84" length
bed width	35"	35"	35"
bed width, with side rails	40.0"	40.0"	40.0"
bed width, with assist bar	36.5"	36.5"	36.5"
overall bed length with head and foot boards	80.5"	84.5"	88.5"
deck height range	7.5" to 28"	7.5" to 28"	7.5" to 28"
safe working load	600 lbs	600 lbs	600 lbs

NOA Twin Elite Riser

The NOA Twin Elite Riser™ Bed incorporates an extra-wide sleep surface for a new level of resident comfort and relaxation that makes them feel like they're back in their own bed. No more tossing and turning in standard-size long-term care beds.

Wide Sleep Deck

All NOA Twin Elite Riser beds come with an extra-wide 39" sleep surface.

Touch a Button

And watch the wide NOA Twin Elite Riser gently rise from 7.5" to 28" for resident safety and caretaker convenience.

Touch Another Button

And watch as the Twin Elite Riser adjusts the head and leg sections of the bed to fit your residents needs, and/or easily adjusts to a comfort chair position.

Guarantee a Good Nights Sleep

With 39" of sleep space your residents will have more freedom to move and turn in bed.



Features Include:

- A 39" wide sleep surface for superior comfort
- High-Low function from 7.5" to 28" for resident safety and caregiver convenience
- 600 pound safe working load
- Four quiet DC motors for smooth operation
- Effortless mobility in high or low positions
- Long linen life, excellent mattress aeration and easy cleaning
- Four Year Warranty on Electronics
- Lifetime warranty on bed deck and frame

Optional Features:

- Safety lock for immobilization
- Head and Foot Boards - Light Oak, Dark Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Sub-acute handset
- Assist Bar, Left, Right and Universal
- Floor Pad, 2" thick, x 2' x 6'
- Therapeutic Foam Mattresses
- Footboard Mounted Staff Control
- Full Overbed Trapeze

ITEM

	80" length	84" length
bed width	39.0"	39.0"
bed width, with side rails	44.0"	44.0"
overall bed length with head and foot boards	84.5"	88.5"
deck height range	7.5" to 28"	7.5" to 28"
safe working load	600 lbs.	600 lbs.

NOA Elite Bariatric Bed

Wide Sleep Deck

All NOA Elite Bariatric Beds come with a comfortable 42" wide sleep deck that easily adjusts from 80" to 84" in length.

Soft - Touch Side Rails

Side Rails are now made with a soft touch vinyl coated steel skeleton. The vinyl coating protects the residents from accidental bruising and the steel skeleton provides a rigid hand hold for ingress/egress and repositioning.

Recliner Position

For added comfort, the bed can be put into a Recliner Position with the bed in a Reverse Trendelenberg position and the knee gatch and head deck raised for greater resident comfort.

Trolley System

NOA Elite Bariatric Beds come with a trolley system which allows movement of the bed from one room to the other or upright storage to save space when not in use.

Lock It In Place

The Elite Bariatric Bed can be moved at any height, at any position or locked in place.



Trolley System



Features Include:

- Soft-Touch Side Rails
- Polyethylene Head and Foot Boards
- Transport Trolley System
- T/RT Capable
- Baseboard Bumper
- Entrapment Compliant

Optional Features:

- Magnetic Lockout Handset
- Assist Bar, Left, Right and Universal
- Floor Pad, 2" thick, x 2' x 6'
- Therapeutic Foam Mattresses
- Footboard Mounted Staff Control
- Full Overbed Trapeze

ITEM	80" length	84" length
bed width	42.0"	42.0"
bed width, with side rails	47.0"	47.0"
overall bed length with head and foot boards	84.5"	88.5"
deck height range	9.5" to 30"	9.5" to 30"
safe working load	800 lbs	800 lbs

NOA Elite R600

Self-Leveling

A built-in sensor in the control box communicates with the actuators to keep the bed level at all times.

600 lb. Safe Working Load

4 Year Warranty

on electronics and lifetime warranty on bed and deck frame.

Hospital Grade Casters

allow for easy transport within the facility.

Fully Adjustable

The Elite R600 provides trendelenberg positioning as well as cardiac chair positioning.



Features Include:

- Height Range : 9.5" at low position to 29.5" at high position
- Standard Hospital Grade 4" Locking Casters
- High density, white Polyethylene head and foot boards with handholds
- Rib Deck design for added durability and increased mattress ventilation
- Can be Moved in any position and at any height
- Goes Straight up and Straight Down

Options Available:

- Assist Bar, Left, Right and Universal
- In-rail Pendant
- Lock Out Hand Set
- Floor Pad, 2" thick, x 2' x 6'
- Therapeutic Foam Mattresses
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- High Impact Plastic head and foot boards
- Standard Laminate Colors available

ITEM	76" length	80" length	84" length
bed width	35"	35"	35"
bed width, with side rails	40.0"	40.0"	40.0"
bed width, with assist bar	36.5"	36.5"	36.5"
overall bed length with head and foot boards	80.5"	84.5"	88.5"
deck height range	9.5" to 29.5"	9.5" to 29.5"	9.5" to 29.5"

NOA LTC Light

Able to effectively meet the needs of your clinical staff, the LTC Light provides a safe working load of 450 lbs. and offers a number of quality standard accessories.

NOA Medical's newest high/low bed system quickly rises from a low of 8.5" to a high of 25".

Fully adjustable

The four motors provides your residents with comfort positioning at the head and foot area, as well as height adjustment.

High-Impact head and foot boards with quick release

Included as standard with the NOA LTC Light adding



to the cost effectiveness of this new bed system from NOA Medical.

Bed system package includes the deck and frame, the electronics, mattress stays at corners, one handset and high-impact head and foot boards in one of five standard wood grain inset colors all for one low price.

Features Include:

- Rises from a low of 8.5" to a high of 25".
- Has 4 DC Motors that provide elevation, head and foot adjustment.
- Has a sturdy flat vented deck providing a solid sleep surface and mattress ventilation.
- Is 35" wide by 80" long.
- Locking Casters.
- Can be moved in any position.
- Easy to clean and disinfect
- Options available:
 - Soft touch side rails
 - Assist bar
 - Foot Board / Staff Controls
 - Standard wood grain head and foot boards



ITEM

80" length

bed width	35"
bed width, with side rails	40.0"
bed width, with assist bars	36.5"
sleep deck length	80.0"
overall bed length with head and foot boards	84.5"
deck height range	8.5" to 25"
safe working load	450 lbs

Accessories

Side Rails




Sub Acute Handset



Safety Lock



Assist Bar



Standard Head and Foot Boards
Available in 5 colors




Impact Resistant Head and Foot Boards with Quick Release
Available in 5 colors



Floor Mats




Color Options



Mahogany




Walnut



Light Oak



Dark Oak



Wild Cherry

Accessories (cont.)



Over Bed Table

The NOA Overbed Table offers solid, durable construction with a progressive design. The adjustment lever provides ease-of-use in raising and lowering the table top, straight up and down from 29" to 43". The sturdy U-Base comes in either black or beige with the tops featuring high pressure laminate available in our standard colors of Light Oak, Dark Oak, Walnut, Wild Cherry and Mahogany. Custom (special order) color choices are also available.



Deck Expansion System

NOA Medical's NEW Deck Expansion System is like no other on the market. Fitting on the NOA Elite Family of Beds, the Deck Expansion System allows you to grow a standard 35" deck to 39" and 42", our twin size decks from 39" to 42" and our bariatric bed system from 42" to 48". Our patent pending design allows for easy attachment to your NOA Bed deck and offers you a quick and easy adjustment for your bed width requirements.

Trapeze

The NOA Trapeze is designed to help patients shift positions while in bed with minimal help.

- Provides safety and security
- Mounts directly onto bed
- Easy to assemble
- Economical
- Fits both standard and bariatric beds
- Rated at 800 lb



The logo for NDA Medical Industries Inc. features the letters 'NDA' in a bold, blue, sans-serif font. The letter 'O' is replaced by a stylized graphic of three vertical white bars of varying heights, creating a circular shape.

MEDICAL INDUSTRIES INC.

801 Terry Lane
Washington, MO 63090
Toll Free: 1-800-633-6068
www.noamedical.com

A metallic, three-dimensional badge with a silver border and a dark blue center. The text 'FINANCING AVAILABLE' is written in white, bold, uppercase letters across the center.

FINANCING AVAILABLE

*Ask us about our
financing program*