



## Elite EX™ Part of the Elite Series of Beds

NOA Medical's Elite EX leads the long term care industry in providing a wider, more comfortable sleep deck for residents.

### Standard Features

- Integrated Deck Expansion System, tool-less expansion of bed from 36" to 39" to 42" width through a simple pull of a button, and 80" to 84" in length.
- Travels from a low of 7" to a convenient high height of 30".
- Unique anti-shear head deck, easily keeps resident from sliding down in bed.
- Center Step Safety Lock
- 600 lb. safe working load
- Self-Leveling motors
- Four year warranty on electronics and 15 year warranty on bed deck and frame
- Swivel locking casters and flip down caster guides
- Baseboard bumper
- 12 button handset

### Optional Features

- Head and Footboards- Light Oak, Wild Cherry, Walnut, Mahogany
- Soft Touch Side Rails
- Standard or Folding Assist Bar
- Pressure Redistribution Foam Mattress
- Footboard Mounted Staff Control
- Full Overbed Trapeze
- Fall Mat



## Technical Specifications

Bed length .....	80" Length (expands easily to 84")
Bed width .....	35" (can be expanded to 39" or 42")
Bed width, with side rails .....	40"
Safe working load .....	600 lbs.
Overall bed length with head and foot boards: .....	84.5"
Deck height range: .....	7" - 30"