



A DIVISION OF NOA MEDICAL INDUSTRIES

Platinum Hospital Transport Stretcher



Introducing the NOA Platinum Patient Transport Stretcher

Our hydraulic patient transport Stretcher offering you multiple standard features at an affordable price.

- Powder coated steel frame and ABS base
- Back board is X-Ray translucent and is fitted with X-Ray cassette holder
- Foldable 5-Staff aluminum alloy guardrail
- Equipped with two USA hydraulic pumps
- Five functions: Adjustable back rest, leg rest, trendelenburg/reverse trendelenburg and height adjustment
- Telescoping IV stand
- 5th Wheel assist for directional control during transport
- 8 inch casters with central locking
- Two retractable push handles for controlled movement
- Mattress with moisture resistant cover
- Oxygen cylinder holder/tray for placement in ABS base
- Rubber Corner bumpers



801 Terry Lane, Washington, MO 63090
855.324.3656
noamedical.com



A DIVISION OF NOA MEDICAL INDUSTRIES

Platinum Hospital Transport Stretcher



Technical Specifications

Feature	Dimension
Overall stretcher width	31"
Overall stretcher length	83"
Low position floor to deck	21.5"
High position floor to deck.	34"
Safe working load	485 lbs.
Caster size	8"
Max head elevation.	80 degrees
Max trend and reverse trend	18 degrees



801 Terry Lane, Washington, MO 63090
855.324.3656
noamedical.com