



A DIVISION OF NOA MEDICAL INDUSTRIES

## Platinum ML Hospital Bed



## Platinum ML Hospital Bed

Economically Priced with Top of the Line Features

- Manual CPR Quick Release Head Deck
- Out of Bed Alarm with Nurse Call Button
- Under Bed Lighting
- Battery Backup
- High Impact Quick Release Head and Footboards
- High Impact Side Rails
- Central Locking System with 5" Casters
- 15 Year Warranty on Deck & Frame Welds
- 4 Year Warranty on Electronics
- Drainage Hooks
- Caster Base Covers
- Self Leveling Motors
- 600 lb. Safe Working Load
- Foot Board Control
- Pillow Speaker



A DIVISION OF NOA MEDICAL INDUSTRIES

801 Terry Lane, Washington, MO 63090

855.324.3656

noamedical.com



A DIVISION OF NOA MEDICAL INDUSTRIES

### Platinum ML Hospital Bed



## Technical Specifications

### Platinum ML Hospital Bed

Feature	Dimension
Frame width . . . . .	41"
Frame length with head and footboard standard length . . . . .	88"
Sleep deck	
width . . . . .	36"
length . . . . .	80"
Bed clearance . . . . .	5.5"

Feature	Dimension
Low position floor to deck . . . . .	11.5"
High position floor to deck . . . . .	31.75"
Max Trend and Reverse Trend . . . . .	14 degrees
Max head elevation . . . . .	68 degrees
Caster size . . . . .	5"



A DIVISION OF NOA MEDICAL INDUSTRIES

801 Terry Lane, Washington, MO 63090  
855.324.3656  
noamedical.com