



A DIVISION OF NOA MEDICAL INDUSTRIES

Platinum NS Hospital Bed



Platinum NS Hospital Bed

Meets the needs of your patients and Clinical Staff at an Affordable Price

- One Button CPR Quick Release Head Deck
- Out of Bed Alarm with Nurse Call Button
- Under Bed Lighting
- Battery Backup
- High Impact Quick Release Head and Footboards
- High Impact Side Rails
- Central Locking System with 5" Casters
- 15 Year Warranty on Deck & Frame Welds
- 4 Year Warranty on Electronics
- Drainage Hooks
- Caster Base Covers
- Self Leveling Motors
- 600 lb. Safe Working Load
- Nurse Panel
- Pillow Speaker
- Meets 60601-2-52



A DIVISION OF NOA MEDICAL INDUSTRIES

801 Terry Lane, Washington, MO 63090

855.324.3656

noamedical.com



A DIVISION OF NOA MEDICAL INDUSTRIES

Platinum NS Hospital Bed



Technical Specifications

Platinum NS Hospital Bed

Feature	Dimension	Feature	Dimension
Frame width	42"	Low position floor to deck	13"
Frame length with head and footboard		High position floor to deck	33.5"
standard length	88"	Max Trend and Reverse Trend	14 degrees
Sleep deck		Max head elevation	68 degrees
width	36"	Caster size	5"
length	80"		
Bed clearance	5.5"		



A DIVISION OF NOA MEDICAL INDUSTRIES

801 Terry Lane, Washington, MO 63090
855.324.3656
noamedical.com