



## Platinum SC Plus Hospital Bed



# Introducing the NOA Platinum SC Plus Bed.

## Quality Features at an Affordable Price

- Patent Pending in Bed Scale System
- Scale never needs calibration
- Patient Hand Pendant
- Transport Handles for ease of transporting bed throughout hospital
- Digital Nurse Panel Providing you real time readout of head deck angle, knee deck angle, height of bed and angle of T/RT.
- Push Button CPR Release Head Deck gently lowering the patient
- 8 Hospital Grade Casters with Central and Directional Locking
- Bed Exit Alarm with Multiple Zones that ties into the Nurse Call System
- Underbed Lighting
- Optional Side Rail Controls for Clinician Ease of Use
- 600 lb. Safe Working Load



801 Terry Lane, Washington, MO 63090  
855.324.3656  
noamedical.com

Platinum SC Plus Hospital Bed



## Technical Specifications

### NOA Platinum SC Plus Hospital Bed

Feature	Dimension	Feature	Dimension
Frame width . . . . .	39.75"	Low position floor to deck . . . . .	11.5"
Frame length with head and footboard (with bumpers) normal length . . . . .	92"	High position floor to deck . . . . .	31.75"
Sleep deck width . . . . .	36"	Max head elevation . . . . .	68 degrees
normal length . . . . .	80"	Max Trend and Reverse Trend . . . . .	14 degrees
Bed clearance under frame . . . . .	5.75"	Caster size . . . . .	5"
Optional Scale System accuracy . . . . .	+ / - 1 % of patient weight		
weigh in any position . . . . .	yes		
scale repeatability . . . . .	2% of patient weight		