



Platinum SCE Plus Hospital Bed



Introducing the NOA Platinum SCE Plus Bed.

Quality Features at an Affordable Price

- Patent Pending in Bed Scale System
- Never needs calibration
- Patient Hand Pendant
- Transport Handles for ease of transporting bed throughout hospital
- Digital Nurse Panel Providing you real time readout of head deck angle, knee deck angle, height of bed and angle of T/RT.
- Push Button CPR Release Head Deck gently lowering the patient
- 8 Hospital Grade Casters with Electronic One Button Push and Central & Directional Locking
- Bed Exit Alarm with Multiple Zones that ties into the Nurse Call System
- Underbed Lighting
- Optional Side Rail Controls for Clinician
- Ease of Use
- 600 lb. Safe Working Load



801 Terry Lane, Washington, MO 63090
855.324.3656
noamedical.com

Platinum SCE Plus Hospital Bed



Technical Specifications

NOA Platinum SCE Plus Hospital Bed

Feature	Dimension	Feature	Dimension
Frame width	39.75"	Low position floor to deck	11.5"
Frame length with head and footboard (with bumpers)		High position floor to deck	31.75"
normal length	92"	Max head elevation	68 degrees
Sleep deck		Max Trend and Reverse Trend	14 degrees
width	36"	Caster size	5"
normal length	80"		
Bed clearance under frame	5.75"		
Optional Scale System			
accuracy	+ / - 1 % of patient weight		
weigh in any position	yes		
scale repeatability	2% of patient weight		