



MEDICAL INDUSTRIES
Hoffmann Family of Companies



NOA Platinum NS Hospital Bed with Prison Modifications

Meeting the needs of your patients and Clinical Staff at an Affordable Price

Features:

- High impact quick release head and footboards
- One hand release, dropdown side rails
- Central and directional locking system with 5" Casters
- 600 LB. safe working load
- 15 year warranty on deck and frame welds
- 4 year warranty on electronics
- Meets 60601-2-52
- Self-Leveling Motors

Prison Modified Features:

- Secure Fasteners on all motors/electronics
- Hidden/Secure wiring designed for safety
- Devices with wiring and plug-ins can be removed and stored until needed
- Bolt on Mattress Retainers
- Central and Directional Locking with removeable locking tool
- Secured Head and Foot Side Rails
- Foley Bag Hooks on both sides of bed
- Lock out available on all bed functions
- Removable Power Cord



MEDICAL INDUSTRIES
Hoffmann Family of Companies



Technical Specifications

NOA Platinum NS Hospital Bed with Prison Modifications

Overall Width (Includes Side Rails).....	39.5"	Safe Working Load.....	600 lbs.
Overall Length (Includes Head & Footboards).....	86"	Max & Reverse Trendelenburg.....	14"
Patient Surface.....	34.5" x 80"	Max Back Angle.....	70"
Bed Clearance.....	5.5"	Max Hip-Knee Angle.....	35"
Deck Height Range.....	17" - 32.5"		



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com