



MEDICAL INDUSTRIES
Hoffmann Family of Companies



NOA Platinum SC Plus Bed with Prison Modifications

Quality Features at an Affordable Price

Features:

- High impact quick release head and footboards
- One hand release, dropdown side rails
- Patented in bed scale system
- Scale never needs calibration
- Bed exit alarm
- Central and directional locking system with 5" Casters
- 600 LB. safe working load
- 15 year warranty on deck and frame welds
- 4 year warranty on electronics
- Meets 60601-2-52
- Self-Leveling Motors

Prison Modified Features:

- Integrated Scale and Bed Exit Alarm
- Secure Fasteners on all motors/electronics
- Hidden/Secure wiring designed for safety
- Devices with wiring and plug-ins can be removed and stored until needed
- Bolt on Mattress Retainers
- Central and Directional Locking with removeable locking tool
- Secured Head and Foot Side Rails
- Foley Bag Hooks on both sides of bed
- Lock out available on all bed functions
- Removable Power Cord



MEDICAL INDUSTRIES
Hoffmann Family of Companies



Technical Specifications

NOA Platinum SC Plus Hospital Bed with Prison Modifications

Overall Width (Includes Side Rails).....	39.5"	Max Back Angle.....	70"
Overall Length (Includes Head & Footboards).....	86"	Max Hip-Knee Angle.....	35"
Patient Surface.....	34.5" x 80"	Scale System Accuracy.....	+/- 1% of patient weight
Deck Height Range.....	17" - 32.5"	Weigh in any Position.....	Yes
Safe Working Load.....	600 lbs.	Scale Repeatability.....	2% of patient weight
Max & Reverse Trendelenburg.....	14"		



801 Terry Lane, Washington, MO 63090
800.633.6068 • noamedical.com